

OL-OL3IT48

Overall Bench Dimensions: 48W x 33D x 37H



Open Instrument Bench with lab-grade chemical-resistant phenolic work surface and inset steel frame with lower shelf on heavy-duty total lock casters, with a 1,200 lb uniform load capacity. Overall bench size: 48W x 33D x 37H.

SPECIFICATIONS

OVERALL SIZE 48W x 33D x 37H	WORK SURFACE 1" lab-grade chemical-resistant phenolic	FRAME Inset fully welded steel, epoxy powder coat finish, with wipe-clean laminate shelf
FRAME SIZE (UNDER TOP) 36W x 31D x 30-1/4H	LOAD CAPACITY 1,200 lb uniform load	CASTERS Heavy-duty total lock

COMPONENTS INCLUDED

QTY	MODEL #	DESCRIPTION
1	OL3IT48	Open Instrument Bench 48W x 33D x 37H with Lab-Grade Chemical Resistant Phenolic Work Surface and Inset Steel Frame (36W x 31D x 30-1/4H) with Lower Shelf on Heavy Duty Total Lock Casters, 1,200 lb Uniform Load Capacity

FINISH OPTIONS



Main image shown in **Gray** (standard). **Also available in White** (shown here). Additional frame colors may be available upon request.