

UNIVERSAL INSTRUMENT BENCH 60W


WITH FULL OPEN WELL

OMNILAB

SOLUTIONS

OL-OL3IT60-FSF-WF

Overall Bench Dimensions: 60W x 33D x 37H



Open Instrument Bench with lab-grade chemical-resistant phenolic work surface and flush steel frame with a full open well for equipment or seating on heavy-duty total lock casters, with a 1,200 lb uniform load capacity. Overall bench size: 60W x 33D x 37H.

SPECIFICATIONS

OVERALL SIZE 60W x 33D x 37H	WORK SURFACE 1" lab-grade chemical-resistant phenolic	FRAME Flush fully welded steel, epoxy powder coat finish
OPEN WELL Full open well for equipment or seating	LOAD CAPACITY 1,200 lb uniform load	CASTERS Heavy-duty total lock

COMPONENTS INCLUDED

QTY	MODEL #	DESCRIPTION
1	OL3IT60-FSF-WF	Open Instrument Bench 60W x 33D x 37H with Lab-Grade Chemical Resistant Phenolic Work Surface and Flush Steel Frame (59-3/4W x 31D x 30-1/4H) with Full Open Well for Equipment or Seating on Heavy Duty Total Lock Casters, 1,200 lb Uniform Load Capacity

FINISH OPTIONS



Main image shown in **Gray** (standard). **Also available in White** (shown here). Additional frame colors may be available upon request.