

COMPLETE QTOF BENCH

WITH SINGLE PUMP ENCLOSURE



OL-MS-QTOF-COMP-S

Overall Bench Dimensions: 72W x 33D x 24H



Enclosed QTof Bench with lab-grade chemical-resistant phenolic work surface and inset steel frame with lower shelf on heavy-duty total lock casters. Includes a freestanding noise suppression module for a single vacuum pump with a mobile base, a phenolic top, acoustic foam paneling, push-pull cooling (220 CFM combined rated exhaust), a vented door, three rear hose grommets, and an oil drip pan with shock absorbers; and a two-drawer bank. Also includes a sit-stand PC workstation (monitor arm and keyboard tray) with a vertical CPU/UPS holder, plus two 120V 8-outlet surge suppressors with 12-foot cords. Overall bench size: 72W x 33D x 24H (Pump Enclosure Interior: 18-1/2W x 25-3/4D x 24H).

SPECIFICATIONS

OVERALL SIZE 72W x 33D x 24H	WORK SURFACE 1" lab-grade chemical-resistant phenolic	FRAME Inset fully welded steel, epoxy powder coat finish, with wipe-clean laminate shelf
PUMP ENCLOSURE Acoustic foam lined, 220 CFM push-pull cooling, thermostat controller	CASTERS Heavy-duty total lock	STORAGE AND PANELS Wipe-clean laminate in a color that coordinates with the frame

COMPLETE QTOF BENCH



WITH SINGLE PUMP ENCLOSURE



COMPONENTS INCLUDED

QTY	MODEL #	DESCRIPTION
1	OL3IT72E-QT	Enclosed Instrument Bench, Reduced Height for QToF Systems 72W x 33D x 24H with Lab-Grade Chemical Resistant Phenolic Work Surface and Inset Steel Frame (60W x 31D x 17-1/4H) with 14 Inch Interior Height Clearance and Lower Shelf on Heavy Duty Total Lock Casters, 1,200 lb Uniform Load Capacity
1	OL3DB16-2D	Drawer Bank 15-3/4W x 22D x 14H with Two (2) Full Extension Lab Drawers (Use with Keyboard Tray)
1	OL3PCW-C	Ergotron® Work Surface Mounted Sit-Stand Workstation, Ergonomic Monitor Arm and Keyboard Tray with Independent Height Adjustment, 20 Inch Height Range, Holds Monitors up to 42 Inch, Slide-Out Left or Right Mouse Tray, Wrist Rest, Integrated Cable Management, Clamp or Grommet Mount, Clamps to Surfaces up to 2.5 Inch Thick (Including Work Surface and Frame)
1	OL3CPU	Ergotron® Side-Mounted Vertical CPU/UPS Holder, Adjustable to Fit Units 2 to 8 Inch Wide, 40 lb (18 kg) Capacity, Open Design for Access to Ports and Power Buttons, Steel and Aluminum Construction
2	OL3E	Tripp Lite® Isobar® Surge Suppressor, Eight (8) 120V NEMA 5-15R Outlets in Isolated Filter Banks for EMI/RFI Noise Protection, 3,840 Joule Rating, 12 Amp Resettable Circuit Breaker (1,440 Watt Capacity), UL Listed, Diagnostic LEDs, All-Metal Housing, 12 Ft Cord, Bench Mounted
1	OL3NSF-S-R	Freestanding Noise Suppression Module for Single Vacuum Pump 22W x 28-1/2D x 27H with SONEX® Acoustic Foam Paneling on Interior, High-Airflow Push-Pull Cooling System with Dual Front Supply Fans and Dual Rear Exhaust Fans (220 CFM Combined Rated Exhaust), Three (3) Flexible Hose Grommets in Rear, Includes Digital Display, Programmable Thermostat Controller with Automatic Fan Speed Control and Alarm Warnings, Door Hinged Right
1	OL3NS-TOP22	Add Lab-Grade Chemical Resistant Phenolic Top to Above Single Pump Noise Suppression Module, 22-1/4W x 28-3/4D, 1/8 Inch Overhang on All Sides Protects Cabinet Edges from Chips and Impacts, Creates Work Surface on Enclosure
1	OL3ODP-XL	Oil Drip Pan for Pump with Mounted Shock Absorbers, Large
1	OL3NSMB-S	Mobile Base for OL3NSF-S 22W x 28-1/2D x 4H

PUMP ENCLOSURE INTERIOR: CLEAR DIMENSIONS



Interior clear dimensions: 18-1/2W x 25-3/4D x 24H (20-1/2D at the top, where the exhaust fans sit at an angle). Confirm pump, hoses, cords and service clearance fit before ordering.

FINISH OPTIONS



Main image shown in **Gray** (standard). **Also available in White** (shown here). Additional frame colors may be available upon request.