

BASIC HPLC CART

OL-HPLC-BAS

Overall Cart Dimensions: 36W x 33D x 37H



Open HPLC Cart with lab-grade chemical-resistant phenolic work surface and inset steel frame with lower shelf on heavy-duty total lock casters. Includes a slide-out waste shelf for waste container storage. Overall cart size: 36W x 33D x 37H.

SPECIFICATIONS

OVERALL SIZE 36W x 33D x 37H	WORK SURFACE 1" lab-grade chemical-resistant phenolic	FRAME Inset fully welded steel, epoxy powder coat finish, with wipe-clean laminate shelf
LOAD CAPACITY 1,200 lb uniform load	CASTERS Heavy-duty total lock	WASTE SHELF Slide-out, chemical-resistant laminate surface

BASIC HPLC CART

COMPONENTS INCLUDED

QTY	MODEL #	DESCRIPTION
1	OL3IT36	Open Instrument Bench 36W x 33D x 37H with Lab-Grade Chemical Resistant Phenolic Work Surface and Inset Steel Frame (30W x 31D x 30-1/4H) with Lower Shelf on Heavy Duty Total Lock Casters, 1,200 lb Uniform Load Capacity
1	OL3WSOS	Open Slide-Out Shelf for Waste Container Storage with Chemical Resistant Laminate Surface 16W x 22D x 1-1/2H

FINISH OPTIONS

Main image shown in **Gray** (standard). **Also available in White** (shown here).



Additional frame colors may be available upon request.