


OL-OL3HALC-3633

Overall Cart Dimensions: 36W x 33D x 24H - 38H



Height Adjustable LC Cart with lab-grade chemical-resistant phenolic work surface, a cutout for routing and protecting waste lines, and a waste shelf below on heavy-duty locking casters. Electronic 2-column leg base (24H - 38H, 300 lb evenly distributed capacity) with a programmable paddle controller with height presets and an anti-collision sensor. Overall cart size: 36W x 33D x 24H - 38H.

SPECIFICATIONS

OVERALL SIZE 36W x 33D x 24H - 38H	WORK SURFACE 1" lab-grade chemical-resistant phenolic with waste line cutout	BASE Electric height-adjustable steel base
CAPACITY 300 lb evenly distributed	CONTROLS Programmable paddle controller, height presets, anti-collision sensor	CASTERS Heavy-duty locking

COMPONENTS INCLUDED

QTY	MODEL #	DESCRIPTION
1	OL3HALC-3633	Height Adjustable LC Cart 36W x 33D x 24H - 38H with Lab-Grade Chemical Resistant Phenolic Work Surface, Cutout in Top for Routing and Protecting Waste Lines, Electronic Height Adjustable 2-Column Leg Base with 300 lb Evenly Distributed Weight Capacity, User Programmable Electronic Paddle Controller with Height Presets, Anti-Collision Sensor, and Phenolic Waste Shelf Below on Heavy Duty Locking Casters

FINISH OPTIONS



Main image shown in **Gray** (standard). **Also available in White** (shown here). Additional frame colors may be available upon request.