

DUAL PUMP ENCLOSURE

WITH MOBILE BASE



OL-NS-DM

Overall Dimensions: 36W x 28-1/2D x 31H



Freestanding Noise Suppression Module for dual vacuum pumps with acoustic foam paneling, push-pull cooling (440 CFM combined rated exhaust), in-door venting, double doors, five rear hose grommets, and two oil drip pans with shock absorbers on a mobile base. Includes a digital display with a programmable thermostat controller, automatic fan speed control, and alarm warnings. Also includes a 6-outlet surge suppressor with a 6-foot cord. Overall size: 36W x 28-1/2D x 31H (Pump Enclosure Interior: 32-1/2W x 25-3/4D x 24H).

SPECIFICATIONS

OVERALL SIZE 36W x 28-1/2D x 31H	PUMP ENCLOSURE Acoustic foam lined, 440 CFM push-pull cooling, in-door venting, double doors	INTERIOR CLEAR 32-1/2W x 25-3/4D x 24H
BASE Mobile base on locking casters	CONTROLLER Programmable thermostat, automatic fan speed, alarms	POWER 6-outlet surge suppressor, 6 ft cord

DUAL PUMP ENCLOSURE



WITH MOBILE BASE



COMPONENTS INCLUDED

QTY	MODEL #	DESCRIPTION
1	OL3NSF-D	Freestanding Noise Suppression Module for Dual Vacuum Pumps 36W x 28-1/2D x 27H with SONEX® Acoustic Foam Paneling on Interior, High-Airflow Push-Pull Cooling System with In-Door Venting, Front Supply Fans and Quad Rear Exhaust Fans (440 CFM Combined Rated Exhaust), Five (5) Flexible Hose Grommets in Rear, Includes Digital Display, Programmable Thermostat Controller with Automatic Fan Speed Control and Alarm Warnings, Double Doors
2	OL3ODP-XL	Oil Drip Pan for Pump with Mounted Shock Absorbers, Large
1	OL3NSMB-D	Mobile Base for OL3NSF-D 36W x 28-1/2D x 4H
1	OL3E-NS	Tripp Lite® Surge Suppressor, Six (6) 120V NEMA 5-15R Outlets, EMI/RFI Noise Filtering, 2,420 Joule Rating, 15 Amp Circuit Breaker (1,800 Watt Capacity), ETL Listed to UL 1449, Diagnostic LEDs, 6 Ft Cord, Bench Mounted

PUMP ENCLOSURE INTERIOR: CLEAR DIMENSIONS



Interior clear dimensions: 32-1/2W x 25-3/4D x 24H (18-1/2H at the rear wall, below the angled exhaust fans). Confirm pumps, hoses, cords and service clearance fit before ordering.

FINISH OPTIONS



Main image shown in **Gray** (standard). **Also available in White** (shown here). Additional frame colors may be available upon request.