


BASIC SEQUENCER BENCH

OL-SEQ-BAS

Overall Bench Dimensions: 60W x 33D x 37H



Enclosed Sequencer Bench with lab-grade chemical-resistant phenolic work surface and inset steel frame with lower shelf on leveling casters for stable, fixed placement. Includes a universal barcode scanner holder. Also includes one 120V 8-outlet surge suppressor with a 12-foot cord. Overall bench size: 60W x 33D x 37H.

SPECIFICATIONS

OVERALL SIZE 60W x 33D x 37H	WORK SURFACE 1" lab-grade chemical-resistant phenolic	FRAME Inset fully welded steel, epoxy powder coat finish, with wipe-clean laminate shelf
LOAD CAPACITY 1,200 lb uniform load	CASTERS Footmaster leveling casters for stable, fixed placement	STORAGE AND PANELS Wipe-clean laminate in a color that coordinates with the frame

BASIC SEQUENCER BENCH

COMPONENTS INCLUDED

QTY	MODEL #	DESCRIPTION
1	OL3IT60E	Enclosed Instrument Bench 60W x 33D x 37H with Lab-Grade Chemical Resistant Phenolic Work Surface and Inset Steel Frame (48W x 31D x 30-1/4H) with Lower Shelf on Heavy Duty Total Lock Casters, 1,200 lb Uniform Load Capacity
1	OL3LC	Leveling Caster Upgrade, Set of Four (4) Footmaster® Leveling Casters to Replace Standard Casters, Tool-Free Ratchet-Style Adjustment Lowers Leveling Pad for Stable, Fixed Placement and Retracts for Mobility, 250 kg (550 lb) Load Capacity per Caster
1	OL3BCSH	Bench Mounted Barcode Scanner Holder, Universal
1	OL3E	Tripp Lite® Isobar® Surge Suppressor, Eight (8) 120V NEMA 5-15R Outlets in Isolated Filter Banks for EMI/RFI Noise Protection, 3,840 Joule Rating, 12 Amp Resettable Circuit Breaker (1,440 Watt Capacity), UL Listed, Diagnostic LEDs, All-Metal Housing, 12 Ft Cord, Bench Mounted

FINISH OPTIONS

Main image shown in **Gray** (standard). **Also available in White** (shown here).



Additional frame colors may be available upon request.