

HEIGHT ADJUSTABLE LAB DESK 60W

WITH DRAWER UNIT



OL-DESK-HA60-DB

Overall Desk Dimensions: 60W x 30D x 29H - 48H



Height Adjustable Lab Desk with lab-grade chemical-resistant phenolic work surface on adjustable leveling glides. Electronic 2-column leg base (29H - 48H, 250 lb evenly distributed capacity) with an electronic controller with height presets and free range adjustment and an anti-collision sensor. Includes a 16W mobile drawer unit with three full-extension drawers, total lock casters, and an anti-tip safety wheel. Overall desk size: 60W x 30D x 29H - 48H.

SPECIFICATIONS

OVERALL SIZE 60W x 30D x 29H - 48H	WORK SURFACE 1" lab-grade chemical-resistant phenolic	BASE Electric height-adjustable steel base on leveling glides
CONTROLS Height presets, anti-collision sensor	CAPACITY 250 lb evenly distributed	STORAGE Wipe-clean laminate in a color that coordinates with the frame

HEIGHT ADJUSTABLE LAB DESK 60W

WITH DRAWER UNIT



COMPONENTS INCLUDED

QTY	MODEL #	DESCRIPTION
1	OL3HAD60	Height Adjustable Lab Desk 60W x 30D x 29H - 48H with Lab-Grade Chemical Resistant Phenolic Work Surface, Electronic Height Adjustable 2-Column Leg Base with 250 lb Weight Capacity (Evenly Distributed), Electronic Controller with Height Presets and Free Range Adjustment, and Anti-Collision Sensor on Adjustable Leveling Glides
1	OL3SF-DB16	Flex Mobile Lab Drawer Unit 15-3/4W x 22D x 30H with Three (3) Full Extension Drawers and Fifth Anti-Tip Safety Wheel Under Bottom Drawer on Total Lock Casters

FINISH OPTIONS

Main image shown in **Gray** (standard). **Also available in White** (shown here).



Additional frame colors may be available upon request.



Also available on wheels

Locking, non-marking nylon casters (wheels, set of four) can replace the leveling glides for an additional charge. Ask us to add model OL3CS to your quote. Casters add 2-1/2" to the overall height, for a range of 31-1/2H to 50-1/2H.