

ESSENTIAL LAB TABLE 60W

WITH DRAWER UNIT



OL-TE60-30-DB

Overall Dimensions: 60W x 30D x 36H



Essential Lab Table with lab-grade chemical-resistant phenolic work surface and fully welded steel frame with front and rear apron on adjustable leveling glides. Includes a 16W mobile drawer unit with a matching phenolic top, three full-extension drawers, total lock casters, and an anti-tip safety wheel. Overall size: 60W x 30D x 36H.

SPECIFICATIONS

OVERALL SIZE 60W x 30D x 36H	WORK SURFACE 1" lab-grade chemical-resistant phenolic	FRAME Fully welded steel with front and rear apron, epoxy powder coat finish
STORAGE Wipe-clean laminate mobile units with phenolic tops, in a color that coordinates with the frame	WHEELS Optional locking stem casters, model OL3TE-MOB	BASE Adjustable leveling glides

ESSENTIAL LAB TABLE 60W

WITH DRAWER UNIT



COMPONENTS INCLUDED

QTY	MODEL #	DESCRIPTION
1	OL3TE60-30	Essential Lab Table 60W x 30D x 36H with Lab-Grade Chemical Resistant Phenolic Work Surface, Fully Welded Steel Frame with Front and Rear Apron on Adjustable Leveling Glides
1	OL3SF-DB16	Flex Mobile Lab Drawer Unit 15-3/4W x 22D x 30H with Three (3) Full Extension Drawers and Fifth Anti-Tip Safety Wheel Under Bottom Drawer on Total Lock Casters
1	OL3SF-TOP16	Add Lab-Grade Chemical Resistant Phenolic Top to Above Flex Mobile Storage Unit, 16W x 22-1/4D, 1/8 Inch Overhang on All Sides Protects Cabinet Edges from Chips and Impacts, Creates Work Surface on Mobile Unit

FINISH OPTIONS

Main image shown in **Gray** (standard). **Also available in White** (shown here).



Additional frame colors may be available upon request.



Also available on wheels

Locking stem casters (wheels, set of four) can replace the leveling glides for an additional charge. Ask us to add model OL3TE-MOB to your quote. Frame and casters together maintain the standard 36H work surface height.