

ESSENTIAL LAB TABLE 60W

WITH HALF SHELF



OL-OL3TE60-30-HS

Overall Dimensions: 60W x 30D x 36H



Essential Lab Table with lab-grade chemical-resistant phenolic work surface and fully welded steel frame with front and rear apron and a half shelf below on adjustable leveling glides. Overall size: 60W x 30D x 36H.

SPECIFICATIONS

OVERALL SIZE 60W x 30D x 36H	WORK SURFACE 1" lab-grade chemical-resistant phenolic	FRAME Fully welded steel with front and rear apron, epoxy powder coat finish
HALF SHELF 15D shelf below the work surface	WHEELS Optional locking stem casters, model OL3TE-MOB	BASE Adjustable leveling glides

COMPONENTS INCLUDED

QTY	MODEL #	DESCRIPTION
1	OL3TE60-30-HS	Essential Lab Table 60W x 30D x 36H with Lab-Grade Chemical Resistant Phenolic Work Surface, Fully Welded Steel Frame with Front and Rear Apron and 15D Half Shelf Below on Adjustable Leveling Glides

FINISH OPTIONS



Main image shown in **Gray** (standard). **Also available in White** (shown here). Additional frame colors may be available upon request.



Also available on wheels

Locking stem casters (wheels, set of four) can replace the leveling glides for an additional charge. Ask us to add model OL3TE-MOB to your quote. Frame and casters together maintain the standard 36H work surface height.