

ESSENTIAL LAB TABLE 60W

WITH POWER BAR



OL-OL3TE60-30-P

Overall Dimensions: 60W x 30D x 36H



Essential Lab Table with lab-grade chemical-resistant phenolic work surface and fully welded steel frame with front and rear apron and integrated backsplash for power (power unit sold separately) on adjustable leveling glides. Overall size: 60W x 30D x 36H (40H at backsplash).

SPECIFICATIONS

OVERALL SIZE 60W x 30D x 36H (40H at backsplash)	WORK SURFACE 1" lab-grade chemical-resistant phenolic	FRAME Fully welded steel with front and rear apron, epoxy powder coat finish
POWER Integrated backsplash; power unit sold separately	WHEELS Optional locking stem casters, model OL3TE-MOB	BASE Adjustable leveling glides

COMPONENTS INCLUDED

QTY	MODEL #	DESCRIPTION
1	OL3TE60-30-P	Essential Lab Table 60W x 30D x 36H (40H at Backsplash) with Lab-Grade Chemical Resistant Phenolic Work Surface, Fully Welded Steel Frame with Front and Rear Apron and Integrated Backsplash for Power (Power Unit Sold Separately) on Adjustable Leveling Glides

Recommended power unit (sold separately): OL3E-16, 16-outlet power distribution strip, 48 inch, 15 ft cord. Add it to your quote as a separate line.

FINISH OPTIONS



Main image shown in **Gray** (standard).
Also available in White (shown here).
Additional frame colors may be available upon request.

Also available on wheels
Locking stem casters (wheels, set of four) can replace the leveling glides for an additional charge. Ask us to add model OL3TE-MOB to your quote. Frame and casters together maintain the standard 36H work surface height.