

ESSENTIAL WORKSTATION 72W

WITH DRAWER UNIT AND COMBO STORAGE



OL-TE72-30-UP-DB-CS

Overall Dimensions: 72W x 30D x 78H



Essential Workstation with lab-grade chemical-resistant phenolic work surface and fully welded steel frame with full-length upper shelving, includes two adjustable shelves with full back and seismic lip, integrated backsplash for power (power unit sold separately), open sides, and open back on adjustable leveling glides. Includes a 16W mobile drawer unit with a matching phenolic top, three full-extension drawers, total lock casters, and an anti-tip safety wheel; and a 16W mobile combo storage unit with a matching phenolic top, one full-extension drawer over a single-door cabinet, and total lock casters. Overall size: 72W x 30D x 78H (work surface at 36H).

SPECIFICATIONS

OVERALL SIZE 72W x 30D x 78H (work surface at 36H)	WORK SURFACE 1" lab-grade chemical-resistant phenolic	FRAME One-piece fully welded steel frame, uprights and backsplash included, epoxy powder coat finish
UPPER SHELVING Laminate shelves on steel brackets with aluminum seismic front lip	STORAGE Wipe-clean laminate mobile units with phenolic tops, in a color that coordinates with the frame	BASE Adjustable leveling glides

ESSENTIAL WORKSTATION 72W

WITH DRAWER UNIT AND COMBO STORAGE



COMPONENTS INCLUDED

QTY	MODEL #	DESCRIPTION
1	OL3TE72-30-UP	Essential Workstation 72W x 30D x 78H (Work Surface at 36H) with Lab-Grade Chemical Resistant Phenolic Work Surface, Fully Welded Steel Frame with Full-Length Upper Shelving, Includes Two (2) Adjustable Shelves with Full Back and Seismic Lip, Integrated Backsplash for Power (Power Unit Sold Separately), Open Sides, and Open Back on Adjustable Leveling Glides
1	OL3SF-DB16	Flex Mobile Lab Drawer Unit 15-3/4W x 22D x 30H with Three (3) Full Extension Drawers and Fifth Anti-Tip Safety Wheel Under Bottom Drawer on Total Lock Casters
1	OL3SF-CS16-L	Flex Mobile Lab Combo Storage Unit 15-3/4W x 22D x 30H with One (1) Full Extension Drawer Over Single Door Cabinet Hinged Left on Total Lock Casters
2	OL3SF-TOP16	Add Lab-Grade Chemical Resistant Phenolic Top to Above Flex Mobile Storage Unit, 16W x 22-1/4D, 1/8 Inch Overhang on All Sides Protects Cabinet Edges from Chips and Impacts, Creates Work Surface on Mobile Unit

Recommended power unit (sold separately): OL3E-16, 16-outlet power distribution strip, 48 inch, 15 ft cord. Add it to your quote as a separate line.

FINISH OPTIONS

Main image shown in **Gray** (standard). **Also available in White** (shown here).



Additional frame colors may be available upon request.



Also available on wheels

Locking stem casters (wheels, set of four) can replace the leveling glides for an additional charge. Ask us to add model OL3TE-MOB to your quote. Frame and casters together maintain the standard 36H work surface height.