


CORE WORKBENCH 60W

OL-OL3WC60

Overall Dimensions: 60W x 30D x 36H



Core Workbench with lab-grade chemical-resistant phenolic work surface and fully welded steel frame with storage supports, enclosed sides, and open back on heavy-duty total lock casters. Overall size: 60W x 30D x 36H.

SPECIFICATIONS

OVERALL SIZE 60W x 30D x 36H	WORK SURFACE 1" lab-grade chemical-resistant phenolic	FRAME Fully welded steel with enclosed sides and storage supports, epoxy powder coat finish
STORAGE-READY Supports for bench-mounted drawer and cabinet units	BACK Open back for easy access to utilities	CASTERS Heavy-duty total lock

COMPONENTS INCLUDED

QTY	MODEL #	DESCRIPTION
1	OL3WC60	Core Workbench 60W x 30D x 36H with Lab-Grade Chemical Resistant Phenolic Work Surface, Fully Welded Steel Frame with Storage Supports, Enclosed Sides, and Open Back on Heavy Duty Total Lock Casters

FINISH OPTIONS



Main image shown in **Gray** (standard). **Also available in White** (shown here). Additional frame colors may be available upon request.