

CORE WORKBENCH 72W

WITH DRAWER UNIT AND COMBO STORAGE



OL-WC72-DB-CS

Overall Dimensions: 72W x 30D x 36H



Core Workbench with lab-grade chemical-resistant phenolic work surface and fully welded steel frame with storage supports, enclosed sides, and open back on heavy-duty total lock casters. Includes a 16W bench-mounted drawer unit with three full-extension drawers and a 16W bench-mounted combo storage unit with one full-extension drawer over a single-door cabinet. Overall size: 72W x 30D x 36H.

SPECIFICATIONS

OVERALL SIZE 72W x 30D x 36H	WORK SURFACE 1" lab-grade chemical-resistant phenolic	FRAME Fully welded steel with enclosed sides and storage supports, epoxy powder coat finish
STORAGE Wipe-clean laminate bench-mounted units in a color that coordinates with the frame	BACK Open back for easy access to utilities	CASTERS Heavy-duty total lock

CORE WORKBENCH 72W

WITH DRAWER UNIT AND COMBO STORAGE



COMPONENTS INCLUDED

QTY	MODEL #	DESCRIPTION
1	OL3WC72	Core Workbench 72W x 30D x 36H with Lab-Grade Chemical Resistant Phenolic Work Surface, Fully Welded Steel Frame with Storage Supports, Enclosed Sides, and Open Back on Heavy Duty Total Lock Casters
1	OL3SC-DB16	Core Lab Drawer Unit 15-3/4W x 22D x 26H with Three (3) Full Extension Drawers, Bench Mounted
1	OL3SC-CS16-L	Core Lab Combo Storage Unit 15-3/4W x 22D x 26H with One (1) Full Extension Drawer Over Single Door Cabinet Hinged Left, Bench Mounted

FINISH OPTIONS

Main image shown in **Gray** (standard). **Also available in White** (shown here).



Additional frame colors may be available upon request.