

CORE WORKBENCH 72W

WITH POWER BAR, DRAWER UNIT, AND COMBO STORAGE

OMNILAB

SOLUTIONS

OL-WC72-P-DB-CS

Overall Dimensions: 72W x 30D x 36H



Core Workbench with lab-grade chemical-resistant phenolic work surface and fully welded steel frame with integrated backsplash for power (power unit sold separately), storage supports, enclosed sides, and open back on heavy-duty total lock casters. Includes a 16W bench-mounted drawer unit with three full-extension drawers and a 16W bench-mounted combo storage unit with one full-extension drawer over a single-door cabinet. Overall size: 72W x 30D x 36H (40H at backsplash).

SPECIFICATIONS

OVERALL SIZE 72W x 30D x 36H (40H at backsplash)	WORK SURFACE 1" lab-grade chemical-resistant phenolic	FRAME Fully welded steel with enclosed sides and storage supports, epoxy powder coat finish
POWER Integrated backsplash; power unit sold separately	STORAGE Wipe-clean laminate bench-mounted units in a color that coordinates with the frame	CASTERS Heavy-duty total lock

CORE WORKBENCH 72W

WITH POWER BAR, DRAWER UNIT, AND COMBO STORAGE

OMNILAB

SOLUTIONS

COMPONENTS INCLUDED

QTY	MODEL #	DESCRIPTION
1	OL3WC72-P	Core Workbench 72W x 30D x 36H (40H at Backsplash) with Lab-Grade Chemical Resistant Phenolic Work Surface, Fully Welded Steel Frame with Integrated Backsplash for Power (Power Unit Sold Separately), Storage Supports, Enclosed Sides, and Open Back on Heavy Duty Total Lock Casters
1	OL3SC-DB16	Core Lab Drawer Unit 15-3/4W x 22D x 26H with Three (3) Full Extension Drawers, Bench Mounted
1	OL3SC-CS16-L	Core Lab Combo Storage Unit 15-3/4W x 22D x 26H with One (1) Full Extension Drawer Over Single Door Cabinet Hinged Left, Bench Mounted

Recommended power unit (sold separately): OL3E-16, 16-outlet power distribution strip, 48 inch, 15 ft cord. Add it to your quote as a separate line.

FINISH OPTIONS

Main image shown in **Gray** (standard). **Also available in White** (shown here).



Additional frame colors may be available upon request.