

FLEX WORKBENCH 60W



OL-OL3WF60

Overall Dimensions: 60W x 30D x 36H



Flex Workbench with lab-grade chemical-resistant phenolic work surface and extra-wide fully welded steel frame with open sides and open back on heavy-duty total lock casters. Overall size: 60W x 30D x 36H.

SPECIFICATIONS

OVERALL SIZE 60W x 30D x 36H	WORK SURFACE 1" lab-grade chemical-resistant phenolic	FRAME Extra-wide fully welded steel, open sides and back, epoxy powder coat finish
STORAGE-READY Room below for mobile drawer and cabinet units	SIDES AND BACK Open for easy access from every side	CASTERS Heavy-duty total lock

COMPONENTS INCLUDED

QTY	MODEL #	DESCRIPTION
1	OL3WF60	Flex Workbench 60W x 30D x 36H with Lab-Grade Chemical Resistant Phenolic Work Surface, Extra-Wide Fully Welded Steel Frame with Open Sides and Open Back on Heavy Duty Total Lock Casters

FINISH OPTIONS



Main image shown in **Gray** (standard). **Also available in White** (shown here). Additional frame colors may be available upon request.