

FLEX WORKBENCH 60W

WITH POWER BAR AND DRAWER UNIT



OL-WF60-P-DB

Overall Dimensions: 60W x 30D x 36H



Flex Workbench with lab-grade chemical-resistant phenolic work surface and extra-wide fully welded steel frame with integrated backsplash for power (power unit sold separately), open sides, and open back on heavy-duty total lock casters. Includes a 16W mobile drawer unit with a matching phenolic top, three full-extension drawers, total lock casters, and an anti-tip safety wheel. Overall size: 60W x 30D x 36H (40H at backsplash).

SPECIFICATIONS

OVERALL SIZE 60W x 30D x 36H (40H at backsplash)	WORK SURFACE 1" lab-grade chemical-resistant phenolic	FRAME Extra-wide fully welded steel, open sides and back, epoxy powder coat finish
POWER Integrated backsplash; power unit sold separately	STORAGE Wipe-clean laminate mobile units with phenolic tops, in a color that coordinates with the frame	CASTERS Heavy-duty total lock

FLEX WORKBENCH 60W

WITH POWER BAR AND DRAWER UNIT



COMPONENTS INCLUDED

QTY	MODEL #	DESCRIPTION
1	OL3WF60-P	Flex Workbench 60W x 30D x 36H (40H at Backsplash) with Lab-Grade Chemical Resistant Phenolic Work Surface, Extra-Wide Fully Welded Steel Frame with Integrated Backsplash for Power (Power Unit Sold Separately), Open Sides, and Open Back on Heavy Duty Total Lock Casters
1	OL3SF-DB16	Flex Mobile Lab Drawer Unit 15-3/4W x 22D x 30H with Three (3) Full Extension Drawers and Fifth Anti-Tip Safety Wheel Under Bottom Drawer on Total Lock Casters
1	OL3SF-TOP16	Add Lab-Grade Chemical Resistant Phenolic Top to Above Flex Mobile Storage Unit, 16W x 22-1/4D, 1/8 Inch Overhang on All Sides Protects Cabinet Edges from Chips and Impacts, Creates Work Surface on Mobile Unit

Recommended power unit (sold separately): OL3E-16, 16-outlet power distribution strip, 48 inch, 15 ft cord. Add it to your quote as a separate line.

FINISH OPTIONS

Main image shown in **Gray** (standard). **Also available in White** (shown here).



Additional frame colors may be available upon request.