

# FLEX WORKBENCH 72W

OL-OL3WF72

Overall Dimensions: 72W x 30D x 36H



Flex Workbench with lab-grade chemical-resistant phenolic work surface and extra-wide fully welded steel frame with open sides and open back on heavy-duty total lock casters. Overall size: 72W x 30D x 36H.

## SPECIFICATIONS

|                                                                 |                                                          |                                                                                       |
|-----------------------------------------------------------------|----------------------------------------------------------|---------------------------------------------------------------------------------------|
| OVERALL SIZE<br>72W x 30D x 36H                                 | WORK SURFACE<br>1" lab-grade chemical-resistant phenolic | FRAME<br>Extra-wide fully welded steel, open sides and back, epoxy powder coat finish |
| STORAGE-READY<br>Room below for mobile drawer and cabinet units | SIDES AND BACK<br>Open for easy access from every side   | CASTERS<br>Heavy-duty total lock                                                      |

## COMPONENTS INCLUDED

| QTY | MODEL # | DESCRIPTION                                                                                                                                                                                |
|-----|---------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1   | OL3WF72 | Flex Workbench 72W x 30D x 36H with Lab-Grade Chemical Resistant Phenolic Work Surface, Extra-Wide Fully Welded Steel Frame with Open Sides and Open Back on Heavy Duty Total Lock Casters |

## FINISH OPTIONS



Main image shown in **Gray** (standard). **Also available in White** (shown here). Additional frame colors may be available upon request.