

FLEX WORKBENCH 72W

WITH DRAWER UNIT AND COMBO STORAGE



OL-WF72-DB-CS

Overall Dimensions: 72W x 30D x 36H



Flex Workbench with lab-grade chemical-resistant phenolic work surface and extra-wide fully welded steel frame with open sides and open back on heavy-duty total lock casters. Includes a 16W mobile drawer unit with a matching phenolic top, three full-extension drawers, total lock casters, and an anti-tip safety wheel; and a 16W mobile combo storage unit with a matching phenolic top, one full-extension drawer over a single-door cabinet, and total lock casters. Overall size: 72W x 30D x 36H.

SPECIFICATIONS

OVERALL SIZE 72W x 30D x 36H	WORK SURFACE 1" lab-grade chemical-resistant phenolic	FRAME Extra-wide fully welded steel, open sides and back, epoxy powder coat finish
STORAGE Wipe-clean laminate mobile units with phenolic tops, in a color that coordinates with the frame	SIDES AND BACK Open for easy access from every side	CASTERS Heavy-duty total lock

FLEX WORKBENCH 72W

WITH DRAWER UNIT AND COMBO STORAGE



COMPONENTS INCLUDED

QTY	MODEL #	DESCRIPTION
1	OL3WF72	Flex Workbench 72W x 30D x 36H with Lab-Grade Chemical Resistant Phenolic Work Surface, Extra-Wide Fully Welded Steel Frame with Open Sides and Open Back on Heavy Duty Total Lock Casters
1	OL3SF-DB16	Flex Mobile Lab Drawer Unit 15-3/4W x 22D x 30H with Three (3) Full Extension Drawers and Fifth Anti-Tip Safety Wheel Under Bottom Drawer on Total Lock Casters
1	OL3SF-CS16-L	Flex Mobile Lab Combo Storage Unit 15-3/4W x 22D x 30H with One (1) Full Extension Drawer Over Single Door Cabinet Hinged Left on Total Lock Casters
2	OL3SF-TOP16	Add Lab-Grade Chemical Resistant Phenolic Top to Above Flex Mobile Storage Unit, 16W x 22-1/4D, 1/8 Inch Overhang on All Sides Protects Cabinet Edges from Chips and Impacts, Creates Work Surface on Mobile Unit

FINISH OPTIONS

Main image shown in **Gray** (standard). **Also available in White** (shown here).



Additional frame colors may be available upon request.