

FLEX WORKSTATION 72W

WITH SPLIT SHELVING

OL-OL3WF72-USP

Overall Dimensions: 72W x 30D x 78H



Flex Workstation with lab-grade chemical-resistant phenolic work surface and extra-wide fully welded steel frame with split upper shelving, includes four adjustable shelves with full back and seismic lip, integrated backsplash for power (power unit sold separately), open sides, and open back on heavy-duty total lock casters. Overall size: 72W x 30D x 78H (work surface at 36H).

SPECIFICATIONS

OVERALL SIZE 72W x 30D x 78H (work surface at 36H)	WORK SURFACE 1" lab-grade chemical-resistant phenolic	FRAME Extra-wide one-piece fully welded steel frame, uprights and backsplash included, epoxy powder coat finish
UPPER SHELVING Split shelving: four laminate shelves on steel brackets with aluminum seismic front lip	POWER Integrated backsplash; power unit sold separately	CASTERS Heavy-duty total lock

FLEX WORKSTATION 72W

WITH SPLIT SHELVING



COMPONENTS INCLUDED

QTY	MODEL #	DESCRIPTION
1	OL3WF72-USP	Flex Workstation 72W x 30D x 78H (Work Surface at 36H) with Lab-Grade Chemical Resistant Phenolic Work Surface, Extra-Wide Fully Welded Steel Frame with Split Upper Shelving, Includes Four (4) Adjustable Shelves with Full Back and Seismic Lip, Integrated Backsplash for Power (Power Unit Sold Separately), Open Sides, and Open Back on Heavy Duty Total Lock Casters

Recommended power unit (sold separately): OL3E-16, 16-outlet power distribution strip, 48 inch, 15 ft cord. Add it to your quote as a separate line.

FINISH OPTIONS

Main image shown in **Gray** (standard). **Also available in White** (shown here).



Additional frame colors may be available upon request.