

FLEX WORKBENCH 72W

WITH POWER BAR



OL-OL3WF72-P

Overall Dimensions: 72W x 30D x 36H



Flex Workbench with lab-grade chemical-resistant phenolic work surface and extra-wide fully welded steel frame with integrated backsplash for power (power unit sold separately), open sides, and open back on heavy-duty total lock casters. Overall size: 72W x 30D x 36H (40H at backsplash).

SPECIFICATIONS

OVERALL SIZE 72W x 30D x 36H (40H at backsplash)	WORK SURFACE 1" lab-grade chemical-resistant phenolic	FRAME Extra-wide fully welded steel, open sides and back, epoxy powder coat finish
POWER Integrated backsplash; power unit sold separately	SIDES AND BACK Open for easy access from every side	CASTERS Heavy-duty total lock

COMPONENTS INCLUDED

QTY	MODEL #	DESCRIPTION
1	OL3WF72-P	Flex Workbench 72W x 30D x 36H (40H at Backsplash) with Lab-Grade Chemical Resistant Phenolic Work Surface, Extra-Wide Fully Welded Steel Frame with Integrated Backsplash for Power (Power Unit Sold Separately), Open Sides, and Open Back on Heavy Duty Total Lock Casters

Recommended power unit (sold separately): OL3E-16, 16-outlet power distribution strip, 48 inch, 15 ft cord. Add it to your quote as a separate line.

FINISH OPTIONS



Main image shown in **Gray** (standard). **Also available in White** (shown here). Additional frame colors may be available upon request.